


BLADELESS TROCAR



BLT 312



One Touch Release Button

- To separate or combine Trocar with one hand
- Easy to separate the sleeve from the main body of trocar assembly
- Easy to take out the tissue



Sleeve Fixation

- Stability and convenience for operation by virtue of a fixation to abdominal wall
- A variety of lengths and shapes in sleeve



Scope Insertion & Bladeless Tip

- It enables to mount the scope on the knife.
- The operator can observe damage to abdominal wall during insertion.
- High safety compared with steel-made blade type
- Precise bladeless tip
- No need heavy strength

VAXCON

BLADELESS
TROCAR

Type

Model

Inner Diameter

Sleeve Length

Sleeve Type

Laparoscopy

Thoracoscopy

BLT 305

● ● ●

5 mm

100mm

BLT 310

● ● ●

10mm

100mm

BLT 312

● ● ●

12mm

100mm

BLT 305 -1

● ● ●

5 mm

75mm

BLT 310 -1

● ● ●

10mm

75mm

BLT 312 -1

● ● ●

12mm

75mm

BLT 405

● ● ●

5 mm

56mm

BLT 410

● ● ●

10mm

56mm

BLT 412

● ● ●

12mm

56mm

BLT 405 -1

● ● ●

5 mm

75mm

BLT 410 -1

● ● ●

10mm

75mm

BLT 412 -1

● ● ●

12mm

75mm

BLT 605

5 mm

56mm

BLT 610

10mm

56mm

BLT 612

12mm

56mm

S

1Cannular

D

1Knife + 2Cannular

T

3Cannular