



Single Patient Use Linear Stapler and Reload

Available of Yellow Reload and 45mm Hand Piece

- White, blue, yellow & green reloads to accommodate to various tissue thicknesses
- 45mm linear stapler, ideal for double stapling technique in low rectal sphincter-preserving procedure



Major Features

- Stapling and closure of intestinal, pulmonary, gastric and vascular incision
- Available of 45 and 60mm models with blue, yellow and green reloads;
30mm model with blue and green reloads
- Unique 30mm vessel stapler with white reload of 3 staple rows, specialized for the best vessel stapling
- One single hand piece can be fired 8 times

Golden Stapler Surgical Co., Ltd.



*Advancing Surgical Technology
to Benefit More People*

Positioning Pin

- Avoid slippage of clamped tissue
- Assist parallel closure
- Ensure complete suturing of the target tissue



Parallel Closure

- Even compression on target tissue
- Provide uniform formation height of all staples
- Effective suturing & hemostasis



Intermediate Position

- Adjust suturing position on target tissue
- Ensure suturing on target position



Separate Closure & Firing Handles

- Clear surgical procedures
- Safe & easy to use



Unique 30mm Vessel Stapler

- White reload, 3 rows of staples, exclusive for vessel stapling
- Best stapling and hemostasis effect



Blue, Yellow, Green & White Reloads

- Accommodate to various tissue thicknesses
- Effective stapling and hemostasis



3 Different Sizes 30, 45 & 60mm

- Accommodate various GI tract widths
- Economic and easy to use



Tactile & Audible Responses

- Safety mechanism for empty reload to prevent misfiring



Code	Colour	Staple Rows	Staple Qty	Closed Staple Height (MM)	Stapling Length (MM)	Applicable Tissue Thickness	Reload Code
LSC30W	White	3	23	1.0	30	Blood vessel, thin tissue	LSCR30W
LSC30B	Blue	2	11	1.5	30	Normal thick tissue	LSCR30B
LSC30G	Green	2	11	2.0	30	Thick tissue	LSCR30G
LSC45B	Blue	2	15	1.5	45	Normal thick tissue	LSCR45B
LSC45Y	Yellow	2	15	1.8	45	Medium thick tissue	LSCR45Y
LSC45G	Green	2	15	2.0	45	Thick tissue	LSCR45G
LSC60B	Blue	2	21	1.5	60	Normal thick tissue	LSCR60B
LSC60Y	Yellow	2	21	1.8	60	Medium thick tissue	LSCR60Y
LSC60G	Green	2	21	2.0	60	Thick tissue	LSCR60G

* Staple from titanium material

Golden Stapler Surgical Co., Ltd.

Add: Building 7A, Jiangsu Wujin Sci-Tech Pioneer Park, 256 Mid Mingxin Road,
Hutang Town, Wujin District, Changzhou, Jiangsu Province, China 213164
Tel: +86 519 85070251 Fax: +86 519 85070252
Web: <http://www.goldenstapler.com>